

Condalab®

Peptones for bioprocess.

Better performance in your culture media.

Application notes



Condalab

Inspired by knowledge

Condalo^w

CHO-S cell assays

- Growth and viability with Condalo^w Soy 3
- Comparison Condalo^w Soy vs Ultrafiltered Peptones 4
- Comparison Condalo^w Meat vs Ultrafiltered/non-ultrafiltered Peptones 5
- Comparison Condalo^w Casein vs Casein peptones 6
- Fed Batch Condalo^w Soy peptone Cell Culture: CHO-S growth and viability 7
- Fed Batch Condalo^w Wheat peptone Cell Culture: CHO-S growth and viability 8
- Yield increase in protein production with Condalo^w Soy 9

HEK-293F cell assays

- Condalo^w vs Ultrafiltered Soy Peptone: HEK-293F Cell Culture Growth 10
- Fed-Batch Growth with Condalo^w Soy 11
- Fed-Batch Growth with Condalo^w Wheat 12

Plasmid production

- Fermentation processes with Condalo^w Soy 13

Condalow® Soy peptone Cell Culture: CHO-S growth and viability.

Abstract

The purpose of this assay is to demonstrate Condalow® Soy peptone is a good animal free supplement for growth and recombinant protein productivity enhancement. In this experiment we focus on cell growth and cell viability.

Introduction

Condalow® Soy Peptone is a product obtained from soy protein which through a controlled manufacturing process ensures a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processing, tissue culture media, vaccines, antibody production, and a wide range of BioPharma processes.


Material and methods

The peptones were prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study we used a CHO-S cell line in a serum-free, chemically defined medium. We used a sterile Thomson 24-well plate; 10.4 ml square well, round bottom, individually wrapped.


We cultured 100,000 cells/ml in 4 ml in duplicates per trial. 1 g/l peptone was added to the CDM and the cells. No peptone was added to the control. The cells were grown for five days in a CO₂ (8%) incubator at 350 rpm and 37°C. Cell densities were determined daily with trypan blue in a Neubauer chamber for 5 days.

Results

Effect of Condalow® Soy Peptone on CHO-S Cell Growth*



Batch to Batch Comparison*



*Chemically Defined Medium

The data presented here are the typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Soy Peptone stimulates the growth of the CHO cells, doubling the number of cells with respect to the control during the first 5 days. Minor variations in cell growth were observed between experimental trials.

Condalow® vs Ultrafiltered Soy Peptone: CHO-S Cell Culture Growth and Viability

Abstract

The purpose of this study is to demonstrate Condalow® Soy peptone is a good animal free supplement for growth and recombinant protein productivity enhancement. In this experiment we focus on cell growth and cell viability and how compare Condalow® Soy Peptone with ultrafiltered (UF) products.

Introduction

Condalow® Soy Peptone is a product obtained from soy protein which controlled manufacturing process ensures a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation process, tissue culture media, vaccines, and antibodies production and a wide BioPharma process.


Material and methods

The peptones were prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study we used a CHO-S cell line in a serum-free, chemically defined medium. We used a sterile Thomson 24-well plate; 10.4 ml square well, round bottom, individually wrapped.

We cultured 100,000 cells/ml in 4 ml in duplicates per trial. 1 g/l peptone was added to the CDM and the cells. No peptone was added to the control. The cells were grown for five days in a CO2 (8%) incubator at 350 rpm and 37°C. Cell densities were determined daily with trypan blue in a Neubauer chamber for 5 days.

Results

Effect of Condalow® and UF Soy Peptone on CHO-S Cell Growth



*Chemically Defined Medium.

The data presented here are the typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Soy Peptone stimulates the growth of the CHO cells, doubling the number of cells compared to the control during the first 5 days. Overall Condalow® Soy Peptone proves to be superior to the ultrafiltered peptones used in this study.

Condalow® vs Ultrafiltered/Non-Ultrafiltered Meat Peptone: CHO-S Cell

Abstract

The purpose of this study is to demonstrate that Condalow® Meat Peptone provides better growth and recombinant protein productivity than Ultrafiltered and non-Ultrafiltered Meat hydrolysates from multiple suppliers. In this experiment we focus on cell growth and cell viability and compare Condalow® Meat Peptone against ultrafiltered (UF) and non-ultrafiltered (NUF) products.

Introduction

Condalow® Meat Peptone is a product obtained from meat protein by a controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibodies production and a wide variety of BioPharma processes.

Material and methods

The peptones were prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study we used a CHO-S cell line in a serum-free, chemically defined medium.


We used a sterile Thomson 24-well plate; 10.4 ml square well, round bottom, individually wrapped.

We cultured 100,000 cells/ml in 4 ml duplicates per trial. 1 g/l peptone was added to the CDM culture. No peptone was added to the control. The cells were grown for six days in a CO₂ (8%) incubator at 350 rpm and 37°C.

Cell densities were determined daily with trypan blue using an automated cell counter.

Results

Effect of Condalow® and UF & NUF Meat Peptone on CHO-S Cell Growth



*Chemically Defined Medium.

The data presented here are the typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Meat Peptone stimulates the growth of the CHO cells, doubling the number of cells with respect to the control during the first 6 days. Overall Condalow® Meat Peptone proves to be a superior nutrient source to ultrafiltered and non-ultrafiltered meat peptones while retaining the low endotoxin benefits of ultrafiltration.

Condalow® Casein Peptone vs Competitors Casein Peptone: CHO-S Cell.

Abstract

The purpose of this study is to demonstrate that Condalow® Casein Peptone provides better growth and recombinant protein productivity than Casein hydrolysates from multiple suppliers. In this experiment we focus on cell growth and cell viability and compare Condalow® Casein Peptone against other Casein Peptone products.

Introduction

Condalow® Casein Peptone is a product obtained from dairy protein by a controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibody production and a wide variety of BioPharma processes.

Material and methods


The peptones were prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study we used a CHO-S cell line in a serum-free, chemically defined medium.

We used a sterile Thomson 24-well plate; 10.4 ml square well, round bottom, individually wrapped.

We cultured 100,000 cells/ml in 4 ml duplicates per trial. 1 g/l peptone was added to the CDM culture. No peptone was added to the control. The cells were grown for six days in a CO₂ (8%) incubator at 350 rpm and 37°C.

Cell densities were determined daily with trypan blue using an automated cell counter.

Results



*Chemically Defined Medium.

The data presented here are the typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Casein Peptone stimulates the growth of the CHO cells, doubling the number of cells with respect to the control during the first 6 days. Overall Condalow® Casein Peptone proves to be a superior nutrient source compared to standard peptones with exceptionally low endotoxin concentration.

Fed Batch Condalow® Soy peptone Cell Culture: CHO-S growth and viability.

Abstract

The purpose of this study is to demonstrate that Condalow® Soy peptone is an effective animal free supplement for growth and recombinant protein productivity enhancement. In this experiment we focus on cell growth and cell viability using our Condalow® Soy peptone as the feed for our cultures.

Introduction


Condalow® Soy Peptone is a product obtained from vegetal protein with controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibody production, and a wide range BioPharma processes.

Material and methods

The peptones were prepared at 20 g/L concentration and sterilized by 0.22 µm filtration. For this case study, we used a CHO-S cell line in a serum-free, chemically defined medium (CDM). We used a sterile 24-well plate; 10.4 mL square well, round bottom.

We cultured 300,000 cells/mL in 4 mL in duplicates per trial. 5 g/L peptones were added as feed to the CDM* at a cell population of 4 million/mL. The cells were grown for five days in a CO₂ (8%) incubator at 350 rpm and 37 C. Viable cell densities were determined daily with trypan blue using an automated cell counter.

Results



*Chemically Defined Medium.

The data presented here are typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Soy peptone stimulates the growth of the CHO-S cells, doubling the number of cells compared with the control at day 4.

Fed Batch Condalow® Wheat peptone Cell Culture: CHO-S growth and viability.

Abstract

The purpose of this study is to demonstrate that Condalow® Wheat peptone is an effective animal free supplement for growth and recombinant protein productivity enhancement. In this experiment we focus on cell growth and cell viability using our Condalow® Wheat peptone as the feed for our cultures.

Introduction


Condalow® Soy Peptone is a product obtained from vegetal protein with controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibody production, and a wide range BioPharma processes.

Material and methods

The peptones were prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study, we used a CHO-S cell line in a serum-free, chemically defined medium (CDM). We used a sterile 24-well plate; 10.4 ml square well, round bottom.

We cultured 300,000 cells/ml in 4 ml in duplicates per trial. 5 g/l peptones were added as feed to the CDM on day 3. The cells were grown for five days in a CO₂ (8%) incubator at 350 rpm and 37°C. Viable cell densities were determined daily with trypan blue using an automated cell counter.

Results



*Chemically Defined Medium.
The data presented here are typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Soy peptone stimulates the growth of the CHO-S cells, doubling the number of cells compared with the control at day 4.

Condalow® Soy Peptone: CHO-S promotion of protein production.

Abstract

The purpose of this study is to demonstrate Condalow® Soy Peptone improves growth and recombinant protein productivity enhancement. In these trials, we focus on protein production.

Introduction

Condalow® Soy Peptone is a product obtained from soy protein with a controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibody production, and a wide range of BioPharma processes.

Material and methods

The peptone was prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this application note, we used a CHO-S cell line in a serum-free, chemically defined medium.

We used a sterile 125ml flask to culture 30ml of cells, we cultured CHO cells in flasks with CDM* (control) and flasks with CDM plus 1g/l Condalow® Soy peptone.

The cells were grown for seven days in a CO2 (8%) incubator at 150 rpm and 37°C.

Cell densities were determined daily with trypan blue using an automated cell counter. Cell viability and population were counted for 4 days post-transfection.

The amount of protein was determined using an FPLC chromatography system.

Results

Condalow® Soy vs Control, cell viability




Fig 1. Viable cell density Condalow® Soy Peptone / Viable cell density Control on the right. Measurements were carried out 4 days post-transfection.

Condalow® Soy vs Control, IgG production




Fig 2. Transfected CHO cells IgG production in presence of Condalow® Soy peptone and control.

Conclusions

In this study we demonstrate that the addition of Condalow® Soy Peptone stimulates the production of IgG protein in CHO-S cells. We observe a 30% increase in protein production with the addition of 1g/l Condalow® Soy Peptone. This data shows an significant increase in antibody yield.

Condalow® vs Ultrafiltered Soy Peptone: HEK-293F Cell Culture Growth.

Abstract

The purpose of this study is to demonstrate the effectiveness of Condalow® Soy Peptone, an animal-free supplement, in enhancing both growth and recombinant protein productivity. The study involves a comparison between Condalow® Soy Peptone and ultrafiltered (UF) products, with a focus on cell growth and viability.

Introduction

Condalow® Soy Peptone is a vegetable protein-derived product with a controlled manufacturing process to ensure a low endotoxin content. This product is a rich source of peptides, vitamins, and carbohydrates, and can be utilized in various applications including fermentation processes, tissue culture media, vaccine production, antibody production, and a wide range of BioPharma processes.


Material and methods

The peptones were prepared at a concentration of 20 g/L and sterilized by 0.22 µm filtration. In this case study, we utilized a HEK-293F cell line in a serum-free, chemically defined medium (CDM). Cells were cultured in a sterile 24-well plate with a 10.4 mL square well, round bottom.

We cultured 600,000 cells/mL in 4 mL duplicates per trial. Peptones were added to the CDM* and the cells at a concentration of 2 g/L. The control group did not receive any peptones. The cells were cultured for eight days in a CO₂ (8%) incubator at 350 rpm and 37°C. Viable cell densities were determined daily with trypan blue using an automated cell counter. The media was completely refreshed on day 4.

Results

Effect of Condalow® and UF Soy Peptone on HEK-293F Cell Growth



*Chemically Defined Medium.

The data presented here are typical values. Some minor variation was observed between assays.

Conclusions

In this study, we demonstrate that the addition of Condalow® Soy Peptone stimulates the growth of HEK-293F cells, resulting in a doubling of the cell count compared to the control by day 7. Our results indicate that Condalow® Soy Peptone provides superior growth compared to the ultrafiltered peptones used in this study.

Fed-Batch Condalow® Soy Peptone: HEK-293F Cell Culture Growth.

Abstract

The purpose of this study is to demonstrate Condalow® Soy peptone is a good animal free supplement for growth and recombinant protein productivity enhancement. In this experiment we focus on cell growth and cell viability.

Introduction

Condalow® Soy Peptone is a product obtained from vegetal protein with a controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibody production, and a wide range of BioPharma processes.

Material and methods

The peptone was prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study we used a HEK-293F cell line in a serum-free, chemically defined medium.

We used a sterile 24-well plate; 10.4 ml square well, round bottom.


We cultured 650,000 cells/ml in 4 ml with duplicates per trial. 2 g/l peptone was added to the CDM and the cells. No peptone was added to the control. The cells were grown for five days in a CO₂ (8%) incubator at 350 rpm and 37°C.

Cell densities were determined daily with trypan blue using an automated cell counter.

The media was completely changed on day 4.

Results

Effect of Condalow® Soy Peptone on HEK-293F Cell Growth*



*Chemically Defined Medium.

The data presented here are the typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Soy Peptone stimulates the growth of the HEK-293F cells, doubling the number of cells compared with the control at day 7. With this feed batch strategy we extend the life of the HEK-293F cells which translates into higher biomass production.

Feed Batch Condalow® Wheat Peptone: HEK-293F Cell Culture Growth.

Abstract

The purpose of this study is to demonstrate Condalow® Wheat peptone is a good animal free supplement for growth and recombinant protein productivity enhancement. In this experiment we focus on cell growth and cell viability.

Introduction

Condalow® Wheat Peptone is a product obtained from vegetal protein with a controlled manufacturing process to ensure a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processes, tissue culture media, vaccines, antibody production, and a wide range of BioPharma processes.

Material and methods

The peptone was prepared at 20 g/l concentration and sterilized by 0.22 µm filtration. For this case study we used a HEK-293F cell line in a serum-free, chemically defined medium.

We used a sterile 24-well plate; 10.4 ml square well, round bottom.


We cultured 650,000 cells/ml with 4 ml in duplicates per trial. 2 g/l peptones were added to the CDM* and the cells. No peptone was added to the control. The cells were grown for five days in a CO₂ (8%) incubator at 350 rpm and 37°C.

Viable cell densities were determined daily with trypan blue using an automated cell counter.

The media was completely changed on day 4.

Results

Effect of Condalow® Wheat Peptone on HEK-293F Cell Growth*



*Chemically Defined Medium.

The data presented here are the typical values. Some minor variation was observed between assays.

Conclusions

In this study we demonstrate that the addition of Condalow® Wheat peptone stimulates the growth of the HEK-293F cells, doubling the number of cells compared with the control at day 7. With this feed batch strategy we extend the life of the HEK-293F cells which translates into higher biomass production.

Condalow® Soy Peptone Fermentation: Plasmid production

Abstract

The purpose of this study is to demonstrate Condalow® Soy peptone is a good animal free substitute in LB fermentation media, matching or increasing plasmid production in a range of competent *E. coli* cells used for cloning.

Introduction

Condalow® Soy Peptone is a product obtained from soy protein which through a controlled manufacturing process ensures a low endotoxin content. This product is an excellent source of peptides, vitamins, and carbohydrates. It can be used in fermentation processing, tissue culture media, vaccines, antibody production, and a wide range of BioPharma processes.

Material and methods

Two competent genetically modified *E. coli* strains were used in this study, JM110 and DH5a.

These cells were transformed with an ampicillin resistance plasmid. The transformant cells were inoculated to 4 ml Casein LB and Soy LB, in duplicates, respectively. The cells were left to grow overnight at 250 rpm at 37°C.

OD was measured the next morning and harvested for miniprep. The plasmid yield was measured with nanodrop.

Results

E. coli Growth and Plasmid production




Fig. 1. In this graph Casein LB OD600 and ng/µl are set to 100% signal to plot both results.

Additional data

DH5a	ng/µl	OD600
Condalow® Soy P	473,711	1.88
Condalow® Soy P	374,134	2.1
Casein P	448,646	1.87
Casein P	438,285	1.85
JM110	ng/µl	OD600
Condalow® Soy P	359,085	1.98
Condalow® Soy P	397,361	1.83
Casein P	234,088	2.06
Casein P	230,739	1.99

Table 1: ng/ul and OD % data, showed in Fig. 1.

Conclusions

For the selected *E. coli* strains the substitution of the Casein peptone for Condalow® Soy Peptone resulted in a very similar OD600 readings and plasmid yield, we can conclude Condalow® Soy peptone is a good substitute for animal derived peptone and in some cases improves the growth and plasmid yield.

Condaflow®

www.condalab.com