

For Research Use

TaKaRa

***E. coli* MV1184 Competent Cells**

Product Manual

I. Components

E. coli MV1184 Competent Cells 100 μ l X 10
pUC119 plasmid (0.1 ng/ μ l) 10 μ l
SOC medium* 1 ml X 10

*: SOC medium	2%	Tryptone
	0.5%	Yeast extract
	10 mM	NaCl
	2.5 mM	KCl
	10 mM	MgSO ₄
	10 mM	MgCl ₂
	20 mM	Glucose

II. Specification

E. coli MV1184 Competent Cells are prepared by Hanahan's method modified by Takara Bio and have a transformation efficiency of $>1 \times 10^7$ cfu/ μ g when 100 μ l of the cells are transformed with 1 ng of pUC119. As *E. coli* MV1184 is an amber suppressor(-) strain in which only non-amber DNA vectors can propagate, this product can be used for selecting amber-mutated DNA. As containing F' plasmid, it can be also used as a host of M13 phage vectors or phagemid vectors in preparation of ssDNA. When transformation of pUC vectors or transduction of M13 phage vector DNAs, recombinants can be selected easily by adding X-Gal and IPTG to a media utilizing α -complementarity to β -galactosidase of the Competent Cells.

X-Gal : 5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside

IPTG : Isopropyl- β -D-thiogalactopyranoside

III. Protocols**A. Transformation with a plasmid vector**

- 1) Thaw *E. coli* MV1184 Competent Cells on ice bath just before use.
- 2) Gently mix cells and transfer 100 μ l of Competent Cells into a 14 ml round-bottom tube.

Important : Do not use a vortex to mix cells.

- 3) Add DNA sample (≤ 10 ng is recommended.)
- 4) Keep in the ice bath for 30 min.
- 5) Incubate cells for 45 sec. at 42°C.
- 6) Return to the ice bath for 1 - 2 min.
- 7) Add SOC medium (pre-incubated at 37°C) up to a final volume of 1 ml.
- 8) Incubate by shaking (160 - 225 rpm) for 1 hour at 37°C.
- 9) Plate on selective media.*
*: 100 μ l or less is recommended for plating on 9 cm diameter dish.
- 10) Incubate overnight at 37°C.

B. Transduction with a M13 phage vector

- 1) Follow step 1) - 8) of III-A.
- 2) Add 200 μ l of the host cells (*E. coli* MV1184, A₆₀₀=0.8 - 1.0) into 3 ml of YT soft agar (pre-incubated at 46 - 48°C).
- 3) Add a proper amount of the solution prepared at step 1) into the agar, mix, and immediately spread it onto a YT-plate.
- 4) Incubate at room temperature for 10 - 15 min. and then, at 37°C overnight.

IV. Precautions for Use

- 1) Place a vial of the Competent Cells in a dry ice / EtOH bath immediately upon removal from -80°C freezer. Keep cells in bath until you are ready to proceed.
- 2) You may use 1.5 ml microcentrifuge tubes instead of 14 ml round-bottom tubes (BD company Code: 352059 or 352057, etc.) for transformation, but it may reduce efficiency.
- 3) When using 100 μ l of Competent Cells, apply 10 ng or less of highly purified plasmid DNA. If not, transformation efficiency might decrease.
- 4) When changing an experiment scale, optimum condition should be considered.
- 5) L-broth or ϕ b-broth can be used instead of SOC medium. In this case, lower efficiency might be obtained.

• <u>L-broth</u> :	<u>Ingredient</u>	<u>per 1 L water</u>
	Bacto tryptone	10 g
	Bacto yeast extract	5 g
	NaCl	5 g
	Adjust to around pH7.5 with 1 N NaOH and autoclave.	

• <u>ϕ b-broth</u> :	<u>Ingredient</u>	<u>per 1 L water</u>
	Bacto tryptone	20 g
	Bacto yeast extract	5 g
	MgSO ₄ • 7H ₂ O	5 g
	Adjust to around pH7.5 with 1 N KOH and autoclave.	

- 6) Medium contents

<u>YT soft agar</u> :	<u>Ingredient</u>	<u>per 100 ml water</u>
	Bacto tryptone	0.8 g
	Bacto yeast extract	0.5 g
	NaCl	0.5 g
	Adjust to around pH7.6 with 1 N NaOH, add agar to the concentration of 0.6%, and autoclave.	

<u>YT-plate</u> :	<u>Ingredient</u>	<u>per liter</u>
	Bacto tryptone	8 g
	Bacto yeast extract	5 g
	NaCl	5 g
	Adjust to pH7.6 with 1 N NaOH, add agar to the concentration of 1.5%, and autoclave.	

- 7) *E. coli* MV1184 cells can be prepared by culturing Competent Cells.
- 8) When adding X-Gal or IPTG, follow the procedures described as below:
 - Add 100 mM IPTG to be 100 μ l/100 ml agar medium and 25 μ l/3 ml soft agar.
 - Add 20 mg X-Gal/ml dimethylhormeamide to be 200 μ l/100 ml medium and 50 μ l/3 ml soft agar.
- 9) It is not recommended to refreeze and store the thawed Competent Cells. However, if necessary, freeze in a dry ice/EtOH bath and return to -80°C . The transformation efficiency can be lowered by more than one magnitude.

V. Transformation efficiency

When 100 μ l of *E. coli* MV1184 Competent Cells were transformed with 1 ng of pUC119 according to III-A and the transformants were selected on a LB-agar plate containing Ampicilin, the transformation efficiency was $> 1.0 \times 10^7$ cfu/ μ g pUC119.

VI. Stability of F' plasmid

The white colonies were less than 1% of the total colonies when *E. coli* MV1184 Competent Cells were transformed with pUC119 and plated on a LB-agar plate containing Ampicillin (100 µg/ml), IPTG (0.3 mM), and X-Gal (60 µg/ml).

VII. Genotype

E. coli MV1184 : *ara*, Δ (*lac-proAB*), *rpsL*, *thi* (ϕ 80 *lacZ* Δ M15), Δ (*srl-recA*) 306:: *Tn10* (*tet^r*)/F' [*traD36*, *proAB⁺*, *lac P*, *lacZ* Δ M15]

VIII. Cell density

1 - 2 X 10⁹ bacteria/ml

IX. Storage

-80°C

Note: If it is not stored at -80°C, transformation efficiency may decrease. In this case, it is recommended to confirm the efficiency by using supplied pUC119 prior to use an application.
Do not store in liquid nitrogen.

X. Reference

1) Hanahan, D. (1983) *J. Mol. Biol.* **166**, 557.

XI. Related product

E. coli MV1184 Electro-Cells (Cat. #9025)
pUC118 DNA (Cat. #3318)
pUC119 DNA (Cat. #3319)
pUC118 DNA/BAP (Cat. #3320 - 3324)
pTV118N DNA (Cat. #3328)
X-Gal (5-Bromo-4-Chloro-3-Indolyl- β -D-Galactoside) (Cat. #9031)
IPTG (Isopropyl- β -D-thiogalactopyranoside) (Cat. #9030)

NOTE: This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc.

Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.
