cloned

Terminal Deoxynucleotidyl Transferase

Code No.2230A

Size: 300 units

Shipping at -20°C

Store at -20°C

Supplied: 5 × TdT Buffer 1 ml

Reagents: 0.1% BSA 1 ml

Lot No.

Concentration: units/ μl

Volume: μl

Expiry Date:

Description :

Terminal deoxynucleotidyl transferase (TdT) does not need a template for its reaction, and catalyzes the incorporation of deoxynucleotides into the 3'-OH termini of single- or double-stranded DNA.¹ It requires an oligodeoxynucleotide of at least three bases as primer.

Storage Buffer:

60 mM	Potassium Phosphate, pH7.2
150 mM	KCl
1 mM	DTT
50%	Glycerol

Source :

E. coli carrying the plasmid which encodes the gene of terminal deoxynucleotidyl transferase

Properties :

- Molecular weight: 60,000
- Optimum pH: 7.2
- Cofactor: This product requires divalent ion of Mg^{2+} , Mn^{2+} or Co^{2+} .
- Substrate specificity

Nucleotide analogs (5-methyl-dCTP, 6-O-methyl-dGTP etc.) can be used as a substrate.

Dideoxythymidine triphosphate and cordycepin triphosphate become terminator.

Unit definition :

One unit is the amount of the enzyme that incorporates 1 nmol of [3H] dATP into acid-insoluble product in 1 hour at 37°C at pH7.2, with calf thymus DNA that is treated with DNase and heat-denatured (activated calf thymus DNA) as initiator.

Reaction mixture for unit definition:

100 mM	Sodium cacodylate, pH7.2
8 mM	MgCl_2
0.1 mM	DTT
0.024%	bovine serum albumin
160 $\mu\text{g/ml}$	activated calf thymus DNA
0.5 mM	[^3H] dATP

Purity :

Nuclease activity is not detected in either of the following cases, as judged from the agarose gel electrophoresis pattern:

1. After incubation of 1 μg of λ DNA-*Hind* III fragments with 15 units of this enzyme for 16 hours at 37°C .
2. After incubation of 1 μg of supercoiled pBR322 DNA with 15 units of this enzyme for 16 hours at 37°C .
3. After incubation of 1 μg of 16s, 23s-rRNA with 15 units of this enzyme for 5 hours at 37°C .

Applications:

- Tailing reactions to add complementary homopolymer tails to vectors and cDNA, such as in the Okayama-Berg method.³
- Labelling the 3'-ends of DNA fragments using dNTP or ddNTP.

Composition of Supplied Reagents (Store at -20°C):

1. 5X TdT Buffer	
500 mM	HEPES, pH7.2
40 mM	MgCl_2
0.5 mM	DTT

2. 0.1% BSA

The supplied buffer has the different composition from that used for unit definition. The supplied buffer has the general composition to be applicable to addition of deoxynucleotide to 3'-OH terminus of DNA.

Application example:

Homopolymeric tailing with dGTP

1. Prepare the following reaction mixture.

DNA	4 - 5 pmol (3'-termini)
5X TdT Buffer	10 μl
0.1% BSA	5 μl
dGTP	0.05-0.5 mM (final)
TdT	10 - 15 U
Sterilized distilled water	
Total	50 μl

2. Incubate the mixture at 37°C for 30 min

3. Add 5 μl of 5 M NaCl and 1 μl of 0.5 M EDTA (pH8.0).

4. Extract DNA with 50 μl of phenol/chloroform/isoamylalcohol (25:24:1) and precipitate DNA with chilled ethanol.

Note :

The reaction is affected by the kinds of bases added (dATP, dTTP, dCTP, dGTP), the kinds of divalent cations in the buffer (Mn^{2+} , Co^{2+} , Mg^{2+}), and the terminal structure of DNA (3'-protruding -OH end, blunt end, 3'-recessed -OH end).²⁾ Thus, it is necessary to find the optimum conditions for each reaction.

The optimal molar ratio of dNTPs and DNA when there are 15-40 nucleotides is 20:1 with 3'-protruding -OH ends, and 100:1 with 3'-recessed -OH ends.

References:

- 1) Bollum, F.J. (1974) in *The Enzymes* (Boyer, P.D., ed.), **10**, 145-171
- 2) Deng, G. and Wu, R. (1983) *Methods in Enzymology*, **100**, 96-116.
- 3) Okayama, H. and Berg, P. (1982) *Mol. Cell. Biol.*, **2**, 161-170.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

cloned

Terminal Deoxynucleotidyl Transferase

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