2X PCR Solution

Premix Ex Taq™ Hot Start Version

Code No. RR030A

Size: 500 μ l x 5

(for 100 PCR reactions)

Shipping at -20°C

Store at -20°C

Storage:

Repeated freeze-thaw cycles may decrease the enzyme activity. Once it thawed, dispense into PCR tubes and store at -20°C . (ex. For 50 μ l PCR reaction, dispense by 25 μ l each tube.)

Description:

This product is an optimized mixture composed of enzyme (*TaKaRa Ex Taq*® HS), reaction buffer and dNTP Mixture as 2-fold concentrations. *TaKaRa Ex Taq* HS is designed to be suitable for Hot Start PCR. It is derived from *TaKaRa Ex Taq* and neutralizing monoclonal antibody to *Taq* DNA polymerase. Non-specific amplification due to mispriming and/or formation of primer dimer before thermal cycling can be prevented, since the antibody inhibits the polymerase activity by binding to the *Taq* DNA polymerase until the temperature elevates. This enzyme can be used in general PCR conditions, since monoclonal antibody is denatured in the initial DNA denaturation step.

Contents :

<i>TaKaRa Ex Taq</i> HS*	: 1.25 units/25 μ l
dNTP Mixture	: 2X conc.; ea. 0.4 mM
<i>Ex Taq</i> Buffer	: 2X conc.; including 4 mM Mg^{2+}

* : Specification of *TaKaRa Ex Taq* HS (Cat. #RR006A)

Unit definition :

One unit is the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

25 mM	TAPS (pH9.3 at 25°C)
50 mM	KCl
2 mM	MgCl_2
0.1 mM	DTT
200 μ M	each dATP·dGTP·dCTP
100 μ M	[^3H]-dTTP
0.25 mg/ml	activated salmon sperm DNA

Purity :

Nicking, endonuclease and exonuclease activity were not detected after the incubation of 0.6 μ g of supercoiled pBR322 DNA, 0.6 μ g of λ DNA or 0.6 μ g of λ -*Hind* III digest with 10 units of this enzyme for 1 hour at 74°C .

Applications:

For DNA amplification by Hot Start PCR.

PCR products:

As most PCR products amplified with *TaKaRa Ex Taq* HS have one A added at 3'-termini, the obtained PCR product can be directly used for cloning into T-vector. Also it is possible to clone the product in blunt-end vectors after blunting and phosphorylation of the end.

PCR test :

- Good performance of DNA amplification by PCR was confirmed by using λ DNA as the template (amplified fragment : 20 kb).
- Good performance of DNA amplification of a β -globin gene by PCR was also confirmed by using human genomic DNA as the template (amplified fragment : 17.5 kb).

General reaction mixture for PCR (total 50 μ l) :

<i>Premix Ex Taq</i> Hot Start Version*	25 μ l
Template	<500 ng
Primer 1	0.2 - 1.0 μ M (final conc.)
Primer 2	0.2 - 1.0 μ M (final conc.)
Sterilized distilled water	up to 50 μ l

* : Please mix gently to be uniform and then use.

PCR conditions :

This enzyme can be used in general PCR conditions, since the monoclonal antibody is denatured in the initial DNA denaturation step. No need for a special step to denature the antibody to *Taq* polymerase.

(Example) Amplification of 1 kb DNA fragment

98 $^{\circ}\text{C}$ 10 sec.] 30 cycles	or	98 $^{\circ}\text{C}$ 10 sec.] 30 cycles
55 $^{\circ}\text{C}$ 30 sec.		68 $^{\circ}\text{C}$ 1 min.		
72 $^{\circ}\text{C}$ 1 min.				

(Note) Denaturation conditions vary depending on the thermal cycler and tubes used for PCR. The recommendation is for 5 - 10 sec. at 98°C , or 20 - 30 sec. at 94°C .

TaKaRa Ex Taq is a registered trademark of TAKARA BIO INC.

Premix Ex Taq is a trademark of TAKARA BIO INC.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC.

If you require licenses for other use, please contact us by phone at +81 77 565 6973 or from our website at www.takara-bio.com.

Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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