

Detection of *Escherichia coli* serogroups DNA using real-time PCR

VALIDATION REPORT

ASSAY

Escherichia coli O157, O26, O111, O103, and O145 detection by 5' nuclease real time-PCR assay.

VALIDATED EQUIPMENT

Applied Biosystems, 7500 model; Thermo Scientific, model PikoReal.

SAMPLES (N)

Enrichments prepared from suitable food samples.

VALIDATION ASSAYS

Start – 01-10-2013

Finish – 01-05-2014

SCOPE OF THE METHOD

The test is designed for the detection the genes associated with the serogroups O157, O111, O26, O103, and O145, based on a 5' nuclease real-time PCR reactions. The method's performance was tested according ISO 22118:2011 and ISO/TS 13136:2012 and fulfils the requirements set in this International Standard.

Conventional serotyping is a method based in antibody-antigen agglutination reactions. These is a very complex and time-consuming method. Recently, rapid and sensitive methods have been found to be attractive alternatives to the *E. coli* conventional serotyping method. These methods are PCR-based methods, in particular real-time-PCR-based, that allows the rapid detection of different antigen genes, which in turn allow the determination of *E. coli* serogroups.

SCIENTIFIC BASIS

The present method is based on a 5' nuclease real-time PCR assay, in which specific fluorescent probes are used to detect the amplified DNA by hybridizing with amplicons. These probes are linked to a fluorophore on one end and a quencher that suppresses fluorescence on the other. If the target sequence is present during the PCR, amplification occurs and the probe is degraded, resulting in fluorescence increase. Fluorescence is measured by a detector and the associated software plots the fluorescence intensity versus number of cycles, allowing the determination of the presence or absence of the target organism.

The method also includes a DNA Internal Control (IC). This control is amplified at the same time as the target DNA sequence, but using a different set of primers and a probe labeled with a second fluorophore. The inclusion of the IC in each reaction avoids false negatives due to the presence of substances inhibitory to PCR and validates any negative result.

SAMPLES

For the validation of the present test were used 5 food samples enriched according the method described in ISO 16654:2001 and 45 *Escherichia coli* strains. After enrichment, STEC colonies suspected were isolated, resuspended in the lysis solution and DNA extracted. These food samples included 5 different matrices, aleatory taken from commercial shops:

- Cheese
- Milk
- Meat
- Processed meat
- Soil

PERFORMANCE CHARACTERISTICS OF THE METHOD

1. ANALYTICAL CONTROLS

Each set tests performed simultaneously with the present method include the following controls:

External positive control (PC): reaction containing mustard DNA

Negative control (NC): reaction without any target DNA

Each individual reaction includes:

Internal amplification control (IC): all reactions include an IC DNA and the corresponding set of primes and probe, labeled with a second fluorophore. This control is amplified at the same time, but independently from target DNA sequence. The inclusion of the IC in each reaction avoids false negatives due to the presence of substances inhibitory to PCR and validates any negative result.

To validate the tests, the controls must show the following results:


| | Target detection Channels 1 | Target detection Channel 2 |
|------------------|--------------------------------|-------------------------------|
| Negative Control | Negative | Positive |
| Positive Control | Positive | Not significant |

If the controls do not match these results, the experiment must be repeated.

The results obtained by the present method shall be interpreted as follows:

| Target detection Channels 1 | Target detection Channels 2 | Interpretation |
|--------------------------------|--------------------------------|----------------|
| Positive | Not significant | Positive |
| Ct = NA | Positive | Negative |
| Ct = NA | Ct = NA | Inhibition** |

**When both target and IC detection are negative, it indicates the presence of an amplification inhibitor and the sample must be tested again after 1/10 dilution.



2. INSTRUMENTS

Real Time instruments must be equipped with FAM and ROX detection channels. No special requirements are necessary of the other instruments required (centrifuge, heating block and micropipettes).

The method was validated in an ABI 7500 instrument (Applied Biosystems) and PikoReal (Thermo Scientific).

3. SPECIFICITY

For this propose we used pre- enriched positive samples for O157, O26, O111, O103, and O145 of *Escherichia coli* that were identified by an accredited method. A total of 5 pre-enriched samples were tested (Table 3.1).

| Code | Strain identification | Other informations |
|-------|--------------------------------|--------------------|
| E 2 | <i>Enriched cheese</i> | |
| E 112 | <i>Enriched milk</i> | |
| E 218 | <i>Enriched processed meat</i> | |
| E 600 | <i>Enriched soil</i> | |
| E 662 | <i>Enriched meat</i> | |

Table 3.1 – List of *Escherichia coli* strains used for inclusivity test

Result:

All strains were detected with the present method, corresponding to 100% inclusivity.

3.2. Exclusivity test

The strains for exclusivity test were chosen following the recommendations of International Standard ISO 22118:2011. A total of 30 non target microorganisms related with the target or occurring in the same habitat and that could cause interference with test results were de tested (Table 3.2). The strains were tested using 10 ng of genomic DNA. The suitability of the DNA extracts for amplification was confirmed amplification of the 16S rRNA bacterial gene, using universal primers. The identification of all strains was confirmed by 16S rRNA gene sequencing. All strains were tested in triplicate.

| Code | Strain identification | Other informations |
|--------|--|---------------------|
| MB 385 | <i>Listeria monocytogenes</i> | CECT 4031 |
| MB 388 | <i>Listeria monocytogenes</i> | Serotype 1/2a ou 3a |
| MB 381 | <i>Listeria innocua</i> | CECT 910 |
| MB 357 | <i>Escherichia coli</i> | NCTC 9001 |
| MB 372 | <i>Campylobacter jejuni subsp jejuni</i> | NCTC 11351 |
| MB 373 | <i>Campylobacter coli</i> | NCTC 11366 |
| MB 399 | <i>Salmonella Typhimurium</i> | CECT 443 |
| MB 22 | <i>Vibrio cholerae</i> | Local isolate |
| MB 356 | <i>Vibrio parahaemolyticus</i> | NCTC 10885 |
| MB 363 | <i>Bacillus cereus</i> | IFM1600 |
| MB 115 | <i>Citrobacter freundii</i> | Local isolate |
| MB 451 | <i>Cronobacter muytjensii</i> | ATCC 51329 |
| MB 141 | <i>Lactobacillus plantarum</i> | CECT 748 T |
| MB 144 | <i>Lactobacillus paracasei subsp paracasei</i> | CECT 4022T |
| MB 150 | <i>Staphylococcus aureus</i> | NCTC 6571 |
| MB 302 | <i>Pseudomonas aeruginosa</i> | Local isolate |
| MB 13 | <i>Serratia marcescens</i> | Local isolate |
| MB 163 | <i>Proteus vulgaris</i> | Local isolate |
| MB 365 | <i>Enterococcus faecalis</i> | Local isolate |
| MB 306 | <i>Enterococcus hirae</i> | Local isolate |
| MB 147 | <i>Staphylococcus epidermidis</i> | Local isolate |
| MB 267 | <i>Shigella flexneri</i> | Local isolate |
| MB 249 | <i>Yersinia enterocolitica</i> | Local isolate |
| MB 178 | <i>Klebsiella pneumoniae</i> | Local isolate |
| MB 449 | <i>Brochothrix thermosphacta</i> | ATCC 11509 |
| MB 274 | <i>Legionella pneumophila</i> | Local isolate |
| MB 278 | <i>Legionella micdadei</i> | Local isolate |
| MB 362 | <i>Clostridium perfringens</i> | Local isolate |
| MF 128 | <i>Saccharomyces cerevisiae</i> | CECT |
| MF 129 | <i>Aspergillus niger</i> | Local isolate |

Table 3.2 – List of strains used for exclusivity test

Result:

The method presented 100% exclusivity as none of the strains showed a positive results and positive result.

4. ROBUSTNESS

Robustness of the present method was determined introducing experimental and technical variations in the parameters described in Table 4.1, using two positive and two negative samples and performed in triplicate.

| Parameter | N° variations | Description |
|-------------------------------|---------------|--------------------------------------|
| Annealing temperature | 2 | + 2°C; - 2°C |
| MgCl2 concentration | 2 | + 0,5mM; - 0,5mM |
| Independent performer | | Not applicable |
| Different Real Time equipment | 2 | ABI 7500 PikoReal, Thermo Scientific |
| Independent laboratories | 1 | Biopremier R&D laboratory |

Table 4.1 – Method's parameters submitted to variations

Result:

Congruent results were obtained for all samples in all parameters variants tested. All samples also presented the expected test results.

5. TRUENESS

Trueness of the method was evaluated using 6 positive and 11 negative food samples. All the samples were tested in triplicate.

Result:

All samples presented the expected test results for all the replicates, corresponding to a trueness of 100%.

Performance declaration

Accordingly with the performed tests and the indicative results obtained, it is considered that the present test allows to achieve suitable results in agreement with the Condalab requirements and is validated.

The present test allows the "Detection of *E. coli* serogroups O157, O26, O111, O103, and O145 using real-time PCR".



Laboratorios Conda S.A.
C/ Forja, 9. Torrejón de Ardoz 28850 Madrid, Spain
T. +34 91 761 02 00