

[Click to view this email in a browser](#)

LifeGuard™

Soil Preservation Solution



Early Access to MO BIO's LifeGuard™ Soil Preservation Solution!

Dear Valued Customer,

MO BIO Laboratories is pleased to offer you a special invitation to our early access program for our groundbreaking new product: LifeGuard™ Soil Preservation Solution (catalog # 12868-100) for stabilization of microbial RNA in soil.

Simply add 4ml LifeGuard Soil Preservation Solution directly to 2g of soil to:

- Achieve high quality gene expression profiles from microbial RNA
- Put your RNA in stasis immediately upon collection
- Maintain soil sample stability for up to 30 days at ambient temperature

Participants will be allowed to purchase 100ml bottles of our new LifeGuard™ Soil Preservation Solution prior to our general release of this product. In return, MO BIO simply requests that you share your experience and observations about the product with our technical team.

[- Act now to protect samples collected during your fall expeditions! -](#)

