

Cryopreserve Beads Blue

Cat. 6082

Maintenance freeze medium for bacterial cultures storage

Practical information

Applications	Categories
Maintenance of bacterial cultures	General use

Industry: General cultivation

Formula in g/L

Beef extract	3	Sodium chloride	5
Tryptone	10	Yeast extract	5
L-Cysteine	0,94	Glycerol (ml)	150

Quality control

Solubility	Appearance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°C)
w/o rests	Liquid	N/A	Pale yellow	7,3±0,2

Microbiological test

Microorganisms	Specification
Escherichia coli ATCC 25922	Good growth
Streptococcus pneumoniae ATCC 49619	Good growth

Storage

Temp. Min.:8 °C
Temp. Max.:25 °C

Bibliography
