

# 1/2 Fraser Broth ISO

Cat. 6703

Enrichment medium for the detection and isolation of *Listeria* in food and environmental samples.

## Practical information

---

Applications	Categories
Selective enrichment	Listeria

Industry: Clinical / Food

Regulations: ISO 11290

## Formula in g/L

---

Enzymatic digest of casein	5	Esculin	1
Beef extract	5	Nalidixic acid	0,01
Potassium dihydrogenphosphate	1,35	Sodium chloride	20
Yeast extract	5	Enzymatic digest of animal tissues	5
Lithium chloride	3	Disodium hydrogen phosphate dihydrate	12
Acriflavine hydrochloride	0,0125		

## Quality control

---

Solubility	Appareance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°C)
------------	------------	--------------------------------	------------------------------	-----------------

## Microbiological test

---

## Storage

---

Temp. Min.:2 °C  
Temp. Max.:25 °C

## Bibliography

---