

WL Nutrient Broth

Cat. 1572

For the determination of microbial flora in beer fermentation processes and manufacturing.

Practical information

Applications	Categories
Selective enumeration	Microorganisms of the brewing industry

Industry: Food / Alcoholic beverages / Fermentation

Principles and uses

WL Nutrient Agar, based on the Green and Grey formulation, is recommended for the control of industrial fermentations, particularly the manufacturing of beer. With a pH of 5,5 true counts of beer yeasts can be made. With a pH of 6,5 the medium is ideal for bakery and distilled spirit yeasts.

Tryptone provides nitrogen, vitamins, minerals and amino acids essential for growth. Dextrose is the fermentable carbohydrate providing carbon and energy. Yeast extract is a source of vitamins, particularly of the B-group. Monopotassium phosphate is the buffer. Potassium, Calcium and Ferric chlorides all provide essential ions for the osmotic balance. Manganese and magnesium sulfate are a source of divalent cations. Bromocresol green is the pH indicator.

Formula in g/L

Dextrose	50	Bromocresol green	0,022
Calcium chloride	0,125	Ferric chloride	0,0025
Magnesium sulfate	0,125	Manganase sulfate	0,0025
Monopotassium phosphate	0,55	Potassium chloride	0,425
Tryptone	5	Yeast extract	4

Preparation

Suspend 60,25 grams of medium in one liter of distilled water. Mix well and dissolve by heating with frequent agitation. Boil for one minute until complete dissolution. Sterilize in autoclave at 121 °C for 15 minutes. Cool to 45-50 °C, mix well and dispense into plates.

Instructions for use

Inoculate and incubate at a temperature of 30 °C and observe after 24-48 hours.

Both WL Nutrient (WLN) and WL Differential (WLD) formulae are used in conjunction, as with WLN Agar the bacteria may not be detected unless the number of yeast cells is very small. 1 plate WLN and 2 plates WLD must be used.

- The WLN Agar plate is incubated aerobically for the total plate count of yeasts.
- One of the WLD Agar plates is incubated aerobically for acetic acid bacteria: Flavobacterium, Proteus, thermophilic bacteria and others.
- The second WLD plate is incubated anaerobically for the investigation of lactic-acid bacteria and species of Pediococcus.

All plates are incubated, in general, at 25 °C in the case of beer, and at 30 °C for bakery and malt alcoholic yeasts. Plates are incubated for 2-10 days and up to 2 weeks, according to the flora present. Counts are made at regular intervals during this period.

Quality control

Solubility	Appareance	Color of the dehydrated medium	Color of the prepared medium	Final pH (25°C)
w/o rests	Fine powder	Clear toasted with a blue tint	Light blue	5,5±0,2

Microbiological test

Incubation conditions: (30 °C / 24-48 h).

Microrganisms

Escherichia coli ATCC 25922

Lactobacillus fermentum ATCC 9338

Saccharomyces cerevisiae ATCC 9763

Specification

Moderate growth

Moderate growth

Good growth

Storage

Temp. Min.:2 °C

Temp. Max.:25 °C

Bibliography

Green, S.R. and P.P. Gray 1950. Paper read at American Society of Brewing Chemist Meeting. Wallerstein Lab. Commun 1 2:43. Green, S.R. and P.P. Gray 1950. A differential procedure applicable to bacteriological investigation in brewing. Wallersteia Lab. Commun 13:357.

MacFaddin J.D. 1985. media for isolation cultivation-identification-maintenance of medical bacteria, vol. 1. p. 854-856 Willians Wilkins, Baltimore, MD.