

Specification

Fluid medium used for determining the oxidative and/or fermentative metabolism of Gram negative bacilli (*Enterobacteriaceae*) according to ISO.

Presentation

20 Tubes
Tube 16 x 113 mm
with: 10 ± 0,2 ml

Packaging Details

1 box with 20 tubes, 16 x113 mm glass tubes, ink
labelled and metallic cap.

Shelf Life

12 months

Storage

8-25°C

Composition

Composition (g/l):

Casein peptone.....	2.00
Sodium chloride.....	5.00
D(+) Glucose.....	10.0
Dipotassium phosphate.....	0.30
Bromothymol Blue.....	0.08
Agar.....	3.00

Description /Technique

Description

Using this medium Hugh and Leifson were able to differentiate Gram negative bacteria into three categories: fermentative, oxidative and inactive.

The current formulation has been adapted to the ISO standard described in the Bibliography.

Technique

Proceed according to the current ISO standard.

The inoculation technique is stab inoculation. One tube is covered with a vaseline® layer in order to induce an anaerobic environment that forces the strain to carry out fermentation.

Fermentative organisms produce a large amount of acid in both the tubes, and this is indicated by the yellow colouration of the Bromothymol Blue indicator. Bacteria that utilise an oxidative metabolic pathway carry out this reaction only in the tube without the vaseline ®. Inactive strains do not use sugars and therefore do not induce any change in either tube.

In some instances a slight blue colouration, probably due to alkalinization by peptone degradation, can occur.

Note: pH decrease may appear over time without affecting its performance.

Quality control

Physical/Chemical control

Color : Green

pH: 6.8 ± 0.2 at 25°C

Microbiological control

Inoculate by stabbing with and without vaseline

Microbiological control according to ISO 11133:2014/ Adm 1:2018.

Aerobiosis. Incubation at 35 ± 2°C. Reading at 24 hours.

Microorganism

Ps. aeruginosa ATCC® 27853, WDCM 00025

Escherichia coli ATCC® 25922, WDCM 00013

Salmonella typhimurium ATCC® 14028, WDCM 00031

Salmonella enterica ATCC® 13076, WDCM 00030

Growth

Good growth- O/F: Pos./Neg.

Good growth- O/F: Pos./Pos. Yellow medium

Good growth- O/F: Pos./Pos. Yellow medium

Good growth- O/F: Pos./Pos. Yellow medium

Sterility Control

Incubation 48 hours at 30-35°C and 48 hours at 20-25°C: NO GROWTH

Check at 7 days after incubation in same conditions

Bibliography

- ATLAS, R.M., L.C. PARKS (1993) Handbook of Microbiological Media. CRC Press, Inc. London.
- COWAN, S.T. (1974) Cowan and Steel's manual for the identification of medical bacteria. 2nd ed. Cambridge University Press. UK.
- DOWNES, F.P. & K. ITO (2001) Compendium of Methods for the Microbiological Examination of Foods. 4th ed. APHA. Washington. DC. USA.
- FORBES, B.A., D.F SAHM & A.S. WEISSFELD (Eds.) (1998) Bailey & Scott's Diagnostic Microbiology. 10th ed. Mosby. St Louis. MO. USA.
- HORWITZ, W. (2000) Official Methods of Analysis of the AOAC International 17th ed. Gaithersburg. MD. USA.
- HUGH, R. & E. LEIFSON (1953) The taxonomic significance of fermentative vs. Oxidative metabolism of carbohydrates by various Gram negative bacteria. J. Bact 66:24.
- ISENBERG, H.D. (1992) Clinical Microbiology Procedures Handbook. ASM Press. Washington. DC. USA.
- ISENBERG, H.D. (1998) Essential Procedures for Clinical Microbiology. ASM Press. Washington. DC. USA.
- ISO 21528-1:2017 Standard. Microbiology of food chain - Horizontal methods for the detection and enumeration of Enterobacteriaceae - Part 1: Detection of Enterobacteriaceae.
- ISO 21528-2 (2017) Microbiology of the food chain – Horizontal methods for the detection and enumeration of Enterobacteriaceae. – Part 2: Colony-count method.
- MURRAY, P.R., E.J. BARON, J.H. JORGENSEN, M.A. PFALLER & R.H. YOLKEN (Eds.) (2003) Manual of Clinical Microbiology. 8th ed. ASM Press. Washington. DC. USA.
- MacFADDIN, J.F. (1985) Media for Isolation-cultivation-identification-maintenance of Medical Bacteria. Vol. I. Williams & Wilkins. Baltimore. MD. USA.
- US FDA (Food and Drug Administrations) (1998) Bacteriological Analytical Manual. 8th ed. AOAC International. Gaithersburg. MD. USA.
- WINN, W., S. ALLEN, W. JANDA, E. KONEMAN, G. PROCOP, G. WOODS & P. SCHRECKENBERER (2006) Koneman's Color Atlas and Textbook of Diagnostic Microbiology. 6th ed. Lippincott, Williams & Wilkins. Philadelphia. PA. USA.